

First Aid Cervical Collars



ONE PIECE FIRST AID COLLAR

Made in polyethylene with washable internal part made in EVA (kind of rubber) fixed to structure by plastic rivets. With a tracheal hole, removable chin support and velcro fastener. X-ray translucent. Available in 6 sizes differentiated by safety closure colour or in a set of 6 in a bag.



TWO PIECES FIRST AID COLLAR

Made in polyethylene with washable internal part made in EVA (kind of rubber) fixed to structure by plastic rivets. With a tracheal hole, removable chin support and velcro fastener on the back. X-ray translucent. Available in 6 sizes differentiated by safety closure colour or in a set of 6 in a bag.

GIMA code	One piece first aid neck collars
34620	One-piece collar - infant
34621	One-piece collar - pediatric
34622	One-piece collar - stout (bull-neck)
34623	One-piece collar - adult small
34624	One-piece collar - adult medium
34625	One-piece collar - adult large - tall
34626	Set of 6 collars in a bag

GIMA code	Two pieces first aid neck collars
34630	Two-piece collar - infant
34631	Two-piece collar - pediatric
34632	Two-piece collar - stout (bull neck)
34633	Two-piece collar - adult small
34634	Two-piece collar - adult medium
34635	Two-piece collar - adult large - tall
34636	Set of 6 collars in a bag

Ambu® Perfect™ ACE



Adjustable in 12 sizes

34640



34641
Adjustable in 16 sizes

AMBU® COLLARS

Ambu® Perfect™ ACE collars are rigid collars, available in 2 versions, adult and pediatric, both adjustable in different sizes/positions.

Adult has 16 adjustments, covering neckless, short, regular and tall. Paediatric has 12 adjustments, covering infant, child, small adult. Made in polyethylene.

GIMA code	Adjustable neck collars
34640	Ambu collar paediatric-yellow
34641	Ambu collar adult-white

Universal Collars



UNIVERSAL CERVICAL COLLARS

Collars made of foam coated PVC. Washable, non toxic, with velcro straps. X-Ray translucent. Available in 3 sizes or in a set of 3 in a bag.

GIMA code	Neck collars
34050	Set of 3 collars in a bag
34051	Cervical collar - small
34052	Cervical collar - medium
34053	Cervical collar - large

Splinter Removal Kit



Small metallic splinters, eye tool device

34190



Anti-abrasive nylon loop

34192



"NEODIMIO" magnetic point

Code 34190 SPLINTER REMOVAL KIT

Includes everything to remove splinters and for first aid medications. In a polypropylene box size 235x190x43 mm:
1 Eye wash solution 100 ml with an eye glass
1 magnifying glass, 2 sterile gauze pads 18 x 40 cm
1 Professional removal pencil, 1 mirror.
1 stainless steel tweezers 11.5 cm

Code 34192 PROFESSIONAL REMOVAL PENCIL

Small metallic splinters eye tool remover. On one side you have a NEODIMIO magnetic point, which enables you to withdraw all metallic splinters avoiding any dangerous eventual eye approach. On the other side you have an antiabrasive nylon loop which is used for non metallic splinters. The body is made of anodised aluminium and the neoprene seals guarantee a better pencil hold. Supplied in a box and CE marked.