

## COLOR



WHITE



SILVER



# A 410

## PISTON COMPRESSOR NEBULIZER

CE 0081 EN 184508

COMPACT DESIGN  
VERY SILENT

### PRODUCT SHEET

APPROVALS •  

SHADED POLE INDUCTION MOTOR • 230V 50Hz 90VA (after voltage on demand)

NEBULISATION RATE • 0.50 ml /min

DRUG RESIDUES • 0.3 ml

MAX AIR PRESSURE • 2.3 + 2.5 bar

MAX AIR FLOW • 13 + 15 l/min

NOISE LEVEL • 48 dB

WEIGHT • 1.36 kg

