

Seca Professional Scales



CE 0108
M

Dial with BMI

seca

27284

Code 27284 SECA 756 - mechanical with height meter

Tradition combined with modern design. With a very sensitive precision mechanism to show weight by the aid of a highly visible red pointer on the circular dial.

User manual in English, Spanish, French, German, Italian.

Technical specifications
Capacity: 160 kg
Sensitivity: 500 g
Weight: 18 kg
Dimensions: 628 x 360 x 1000 mm
Platform: 360 x 350 mm with wheels



CE 0108
M

27296

Code 27296 NEW SECA 711 - mechanical Code 27298 HEIGHT METER


Mechanical column scale with weight adjustment at eye level.

Weight can be measured in a conventional way: by aligning two chrome plated indicators.

User manual in English, Spanish, French, German, Italian.

Technical specifications
Capacity: 220 kg
Sensitivity: 100 g
Weight: 23.5 kg
Dimensions: 520 x 1556 x 520 mm
Platform: 350 x 360 mm

Fitness Scale



Large LCD display digits height 17 mm
 USB port opposite LCD for direct connection to PC for data management
 Integrated printer (printout includes date/time, weight, data analysis and evaluations and recommendations for nutrition and fitness activities)
 Glass surface easy to clean

SOEHNLE
PROFESSIONAL

27286

Code 27286 FITNESS SCALE 7850

Its body analysis function can measure muscle, fat and water content, in addition to weight, as well as pulse rate and the current energy conversion.

The recording takes place by hand contact, saving time enabling the user to stand comfortably on the scale with shoes.

In addition, clothing is taken into consideration in the analysis evaluation.

Technical specifications
Capacity: 180 kg
Sensitivity: 100 g (0.2 lbs)
Kg/lbs selection function
Dimensions: 365x710x1290 mm
Platform: 365x485x45 mm
Power supply: 100-240 V, 50/60 Hz

Veterinary Scale



Code 27257 ELECTRONIC VET SCALE

Robust metal platform with a separate digital display.

Technical specifications
Capacity/Sensitivity: 150 kg /50 g 250 Kg/ 100 g
Dimensions: 950 x 500 x 45 mm
Rechargeable battery pack 15V, 300 mA
Full recharging time: 14 hours
Operating time from full recharge: 8 hours
Power supply: 230V/50Hz

Chair scale - see page 36

Bed Digital Scale

Code 27270 SOEHNLE BED SCALE
 Scale consists of 4 weighing units, display terminal and trolley. BMI, hold and dialysis function. The weighing units have to be placed under the bed castors from the side. The bed will be lifted by foot activity without any effort. Fits most types of beds with wheels up to 150 mm.

Technical specifications
Capacity: 500 kg /1100 lb
Sensitivity: 0-300 kg: 100 g/0.2 lb
300-500 kg: 200 g/0.4 lb
Dimensions: 860 x 480 x 1050 mm - 25 Kg
Power supply: 230V/50Hz Made in Germany



Digital Scale

Code 27288 PEGASO DIGITAL SCALE with height meter (70x190 cm)
 Professional digital scale supplied with integrated height meter. This scale can be operated by mains or batteries.

CE 0476

Technical specifications
Capacity: 200 kg
Sensitivity: 100 g
Dimensions: 535x275x940 mm
Platform: 372x275 mm
Weight: 15 kg
Power supply: 220 V - 50 Hz
Battery accumulator: 6V 5AH
Recharging time: 4-5 hours
Operating time when fully charged: 36 hours

